



VET Diagnostix

BCV-ROTA-Crypto-EK99-CPE Ag Combo Test

Easy operation
Simple to operate and easy to interpret

Fast detection
Get a result in 10 minutes

Reliable principle
Rapid immunochromatographic technology

High Accuracy
 ROTA Ag: 96.17% EK99 Ag: 96.19%
 BCV Ag: 97.14% CPE Ag: 96.53%
 CRYPTO Ag: 97.13%




Disease Introduction

Bovine Coronavirus, Rotavirus, Cryptosporidium, Escherichia coli K99, and Clostridium perfringens are common pathogens that can cause diarrhea in cattle, especially in calves, which can lead to severe consequences, including death.

Intended Use

The BCV-ROTA-Crypto-EK99-CPE Ag Combo Test is a rapid diagnostic tool for detecting Bovine Coronavirus antigen, Bovine Rotavirus antigen, Bovine Cryptosporidium antigen, E.coli K99 and Clostridium perfringens antigen.

Test Procedure and Result Interpretation



- Positive (+):** Both control line C and test line T have wine-red bands, indicating that the sample contains coronavirus, rotavirus, cryptosporidium, E.coli K99 or clostridium perfringens antigen.
- Negative (-):** The control line C has wine-red band and the test line T has no color, which means that the sample does not contain coronavirus, rotavirus, cryptosporidium, E.coli K99 or clostridium perfringens antigen.
- Invalid:** No colored band appears in C zone. No matter if T band appears.

Product information

Product name	Specimen	Specification	Shelf Life	Storage
BCV-ROTA-Crypto-EK99-CPE Ag Combo Test	feces	10T/box	18 months	2-30°C

E-mail: sales@aijinbio.com
 Tel: +86-571-85366657 / +86-15868809018
 Web: www.vet-diagnostix.com

HANGZHOU EVEGEN BIOTECH CO., LTD.

